

Power-Pro® ELITE

True Low Air Loss with Alternating Pressure Mattress System

POWER-PRO® ELITE - MODEL 9500

The Power-Pro® Elite Mattress System combines two clinically effective therapies: True Low Air Loss and Alternating Pressure, in a portable system. An optional pulsating low air loss therapy feature is available. The mattress utilizes 18 independent 9" deep air cells to provide proper pressure redistribution for the prevention and treatment of Stages I-IV pressure ulcers. The new Power-Pro VW Elite bariatric model combines variable width convenience with durability and wound prevention features. Available in three models to fit the most popular bariatric bed frames. The Power-Pro VW Elite fits 35", 39", 42", 48" and 54" bed frames. The Power-Pro® Elite True Low Air Loss Alternating Pressure Mattress System supports patients weighing up to 1000 pounds. Audible alarm system alerts caregiver of power outages.

Vyvox - III® - Urethane coated, multi-stretch, moisture resistant / vapor permeable, low shear, quilted cover, transfers moisture away from the patient's skin. The cover is washable, anti-microbial, fluid-resistant and fire-retardant.

Ventilated low air loss cells have air holes throughout the mattress to help minimize skin maceration by reducing excess moisture under patient. Air cells are constructed of TPU thermal polyurethane material. Dual bladder design prevents patients from bottoming out, while sitting up and will maintain air during power outages.

CPR Valve quickly deflates mattress system for medical emergencies. A simple twist releases air immediately.



Hose connections feature a transport locking device that ensures the mattress remains inflated for patient transport.



Model 9500
Mattress & Pump System



Power Pro® VW
Expandable Mattress System

The blower pump is controlled by a micro-processor which allows the caregiver total control of comfort and chosen therapy. Cycle time ranges from 5 to 95 minutes. The blower pump inflates the mattress within minutes, reducing the waiting time for therapeutic treatment to begin. Optional pulsating low air loss pump Model 9511-PS allows caregivers to provide pulsating low air loss therapy.



 **Blue Chip**
Medical Products, Inc.

"Setting the Standard"

POWER-PRO®

SPECIFICATIONS


MODEL 9500 - POWER-PRO® ELITE MATTRESS & PUMP

Dimensions: 36" x 80" x 9"
Weight: 23 lbs.
Air Cells: 18 - 9" Air cells
 Ventilated bladders with air holes throughout the mattress.

Mattress Cover: Vyvex-III® - Urethane coated, multi-stretch, low shear, moisture resistant / vapor permeable quilted cover. Removable, zippered cover, is anti-microbial and meets California Technical Bulletin #117 for fire-retardancy.



MODEL 9511 - POWER-PRO® ELITE PUMP ONLY

Dimensions: 17L" x 7W" x 10.5H"
Weight: 12 lbs.
Electrical: 120V 60 hz
Listings: 
Cycle Time: 5 - 95 Minutes
Air Flow: 60 - 250 Liters per minute
Current: 5A
Ground Resistance: 0.1 ohms maximum



Model 9511 Pump only

ORDERING INFORMATION

Model 9500	Power-Pro® Elite	650 lbs. Capacity	36" x 80" x 9"
Model 9542BAR	Power-Pro® Elite	800 lbs. Capacity	42" x 80" x 9"
Model 9548BAR	Power-Pro® Elite	1000 lbs. Capacity	48" x 80" x 9"
Model 9511	Power-Pro® Elite pump only		
Model 9511-PS	Power-Pro® Elite pump with pulsating low air loss		
Model 9502	Power-Pro® Mattress only with Zippered Quilted Cover 36"		
Model 9502-42	Power-Pro® Mattress only with Zippered Quilted Cover 42"		
Model 9502-48	Power-Pro® Mattress only with Zippered Quilted Cover 48"		
Model 9500-EXP-1	Power Pro® Variable Width Size: 35"-39"-42", 800 lbs.		
Model 9500-EXP-2	Power Pro® Variable Width Size: 35"-42"-48", 1000 lbs.		
Model 9500-EXP-3	Power Pro® Variable Width Size: 39"-48"-54", 1000 lbs.		

*To order pulsating low air loss feature add suffix "PS" to each of the mattress models listed..
 Power Pro® Elite meets Group 2 Support Surface requirements for Medicare E0277.
 Patient must meet all medical criteria to qualify for reimbursement.

WARRANTY

Pumps: 1 year limited
Mattress: 1 year limited
Cover: 6 months limited

PDAC Approved
 Medicare E0277



7-11 Suffern Place, Suffern NY 10901

Tel: (845) 369-7535 • Fax: (845) 369-7633 • Toll Free: (800) 795-6115
www.bluechipmedical.com

"Setting the Standard"

©Blue Chip Medical Products, Inc. 2009

- Alternating Pressure Mattresses
- Alternating Pressure Pads & Pumps
- Low Air-Loss Mattress Systems
- Lateral Rotation Mattress Systems
- Gel Mattresses
- Gel Mattress Overlays
- Static-Air Mattress Overlays
- Water Mattress Overlays
- Wheelchair Cushions
- Patient Lift Slings
- Side Rail Bumpers
- Cart Covers
- OR Replacement Table Pads
- Replacement Stretcher Pads
- Cubicle Curtains
- Hospital Pillows