


Institutional Products

REPLACEMENT STRETCHER & OPERATING ROOM MATTRESSES



1 Section



2 Section




3 Section

All replacement mattresses are designed to meet the heavy-duty requirements necessary for institutional use. The foam core is high density polyurethane foam that meets California Technical Bulletin #117 for fire-retardancy. A selection of foam and materials are available to meet more stringent fire standards such as California Technical Bulletin #129, Boston Fire Code 1X-11 and Federal Flammability Standard 16 CFR 1633 when specified. Your choice of square, round or cut corners are available. Thickness availability: 2", 3" or 4". Available in 1, 2 or 3 sections with either a zippered or sealed lock-stitch closure. Velcro straps, fabric hinges and GU cut-outs are available.

The mattresses are manufactured in a variety of cover materials: the conductive material is a UL® registered fabric approved by the Underwriters Laboratories. The non-conductive is a durable bonded vinyl leatherette. Waterfall flaps cover zippers for added protection against infiltration of fluids. Both materials are anti-microbial, non-allergenic and adhere to California Technical Bulletin #117 for fire-retardancy.

PRESSURE REDISTRIBUTION STRETCHER & OPERATING ROOM MATTRESSES



Pressure redistribution is essential for extended surgical procedures. The combination of PPF® (Pressurized Polymer Foam), Soft Flo® Gel and Viscotec® Visco elastic foam reduces interface pressures and shearing forces. Vyvex III stretch cover is anti-microbial and fire retardant.

CALL TO ORDER, PLEASE SPECIFY SIZE & FEATURES :

SECTIONS:	1 • 2 • 3
SIZE:	CUSTOM
THICKNESS:	2" • 3" • 4"
COVER:	CONDUCTIVE • NON-CONDUCTIVE
CORNERS, HEAD:	SQUARE • ROUND • CUT-OFF
CORNERS, FOOT:	SQUARE • ROUND • CUT-OFF
CLOSURE:	ZIPPERED • SEALED
FEATURES:	VELCRO STRAPS • FABRIC HINGE • GU CUT-OUT



EQUIPMENT AND CART COVERS

We manufacture covers for all departments within a facility: emergency room, linen, operating room, crash carts. All covers are customized to fit your original equipment. Blue Chip offers many different fabrics, such as heavy-duty vinyls, clear see-through vinyl, 400 denier nylon as well as heavy duty ballistic nylon. Each cover can be fabricated with velcro, zipper or lock closure. Silk screening is available with facility logo.

SPECIFY WHEN ORDERING

COVER: VINYL • CLEAR VINYL • DENIER NYLON • BALLISTIC NYLON
CLOSURE: VELCRO • ZIPPER • LOCK CLOSURE

MODEL	SIZE
BC-CC	CUSTOM



7-11 Suffern Place, Suffern NY 10901

Tel: (845) 369-7535 • Fax: (845) 369-7633 • Toll Free: (800) 795-6115

www.bluechipmedical.com

©Blue Chip Medical Products, Inc. 2009

"Setting the Standard"

Blue Chip Medical Products Inc.

Patient & Hospital Accessories

INSTITUTIONAL PRODUCTS

- Stretcher Pads
- Surgical Bolster Rolls
- Equipment and Cart Covers



PATIENT ROOM PRODUCTS

- Folding Floor Mats
- Side Rail Bumpers
- Wedge Guards®
- No Gap® Side Guards
- Mattress Extenders



HOSPITAL & PATIENT ACCESSORIES


- Lift Slings
- Plastic Equipment Bags
- Hospital Pillows



"Setting the Standard"

Hospital & Patient Room Accessories


DEFENDER® FLOOR MATS Patient safety is a must!
Floor mats protect the patient during an accidental fall from bed.



The inner core is a fire-retardant high compression polyurethane foam covered with anti-microbial, fire-retardant vinyl, meeting California Technical Bulletin #117. All mats are 2" and are available in non-fold, bi-fold or tri-fold design. Available in 3" or 4" height as a special order along with Velcro straps.

MODEL	SIZE	SPECIFY WHEN ORDERING:
BCFM 2460	24" x 60"	COLOR: RED • GREY
BCFM 2472	24" x 72"	THICKNESS: 2" • 3" • 4"
BCFM 3072	30" x 72"	STYLE: BI-FOLD • TRI-FOLD
BCFM 3672	36" x 72"	NON-FOLD
CUSTOM SIZES AVAILABLE		OPTIONS: NON-SKID BOTTOM


DEFENDER® EDGE FLOOR MATS



The Defender® "Edge" low profile floor mat has beveled edges that allows for easy wheelchair access and helps to reduce accidental tripping. The inner core is a high compression polyethylene foam covered with Blue Chip's Chemsafe® anti-microbial, fire-retardant vinyl with a non-skid base. The Defender® "Edge" low profile floor mat can be stored under the bed or folded when not in use. Meets Cal. Technical Bulletin #117 for fire retardancy.


MODEL	SIZE
BCEDG-2466	24"W x 66"L x 1.5"H
BCEDG-3072	30"W x 72"L x 1.5"H

SIDE RAIL BUMPERS



Patients are protected from accidental injury due to banging or slipping through a bed rail. The inner core is a fire-retardant high-compression polyurethane that is covered with an anti-microbial, fire-retardant vinyl cover, meeting California Technical Bulletin #117. Bumpers are sold in sets and attach to rail with Velcro closures.

MODEL	SIZE
BCM 211	11" x 70" x 1"
BCM 211-2	11" x 70" x 2"
BCM 215	15" x 60" x 1"
BCM 215-2	15" x 60" x 2"
BCM 216	11" x 30" x 2" (for Hill Rom Beds)
CUSTOM SIZES AVAILABLE	



WEDGE GUARD®

The Wedge Guard doubles as a side rail gap protector in addition to being a patient positioner. Patients are safe from falls and accidental slippage through any gaps between the mattress and side rail. High density foam and medical grade vinyl cover material meet California Technical Bulletin #117 for fire-retardancy and is anti-microbial. Available in full and half rail sizes

MODEL	SIZE
BWG36 Half Rail	36"
BWG70 Full Rail	70"
CUSTOM SIZES AVAILABLE	



MATTRESS EXTENDERS

This versatile device is used to extend a mattress for the needs of a larger patient. When a sheet is placed over the mattress and extender, a unified surface is formed. The extender is covered with Chemsafe® medical grade vinyl. Cover and foam are anti-microbial and meet California Technical Bulletin #117 for fire-retardancy. We can customize a unit to fit any bed size.

MODEL	SIZE
BCME126	6" x 36" x 6"
CUSTOM SIZES AVAILABLE	

"Setting the Standard"

Hospital & Patient Room Accessories

LIFT SLINGS



When safety and reliability for your patients is your main concern, our patient lift slings are up to the task. Available in Canvas, Nylon and Mesh. Our manufacturing process ensures guaranteed strength and durability for homecare and institutional use. All slings are interchangeable and can be used with other manufacturers lifts. Slings work with "V" bar frame only.



MODEL	DESCRIPTION	MATERIAL
7500	REPLACEMENT sling	CANVAS
7501	REPLACEMENT sling with COMMODE CUT-OUT	CANVAS
7502	REPLACEMENT sling	NYLON
7503	REPLACEMENT sling with COMMODE CUT-OUT	NYLON
7504	REPLACEMENT sling	MESH
7505	REPLACEMENT sling with COMMODE CUT-OUT	MESH
7506	REPLACEMENT chains with "S" hooks	
7508	BARiatric Sling, blue heavy duty, 550lb capacity	MESH

PLASTIC EQUIPMENT BAGS

Protect your equipment by keeping them dust free after cleaning. Ideal for complying with accreditation. All bags are a heavy duty weight, 1.5 mil plus. Ideal for JACHO and CHAPS.



MODEL	DESCRIPTION	SIZE	ROLL QUANTITY
8001-CW	COMMODE, CONCENTRATOR, WHEELCHAIRS & WALKER	30" x 12" x 45"	100 per roll
8002-MT	FOR MATTRESSES	36" x 7" x 90"	100 per roll
8003-SM	FOR SUCTION, NEBULIZERS, BIPAP AND CPAP	21.5" x 25"	100 per roll

HOSPITAL PILLOWS

All pillows are manufactured using a "hypo-allergenic" high lofted fire-retardant fiber fill. We offer four different pillow tickings: vinyl, breathable vinyl, navy nylon and poly/cotton. The nylon and vinyl covers are anti-microbial, and meet California Technical Bulletin #117 for fire-retardancy.



MODEL	SIZE	COVER
2026-PC	20" x 26"	Poly/COTTON
2026-V	20" x 26"	Vinyl
2026-B	20" x 26"	BREATHABLE Vinyl
2026-N	20" x 26"	Nylon

NO GAP® SIDE GUARDS

Protect your patients from accidental slippage through any gaps between the mattress and side rail. High density foam and medical grade vinyl cover material meet California Technical Bulletin #117 for fire retardancy and are anti-microbial. Must be removed when rails are not in use.

MODEL	SIZE
BCNG72	72L" x 2W" x 8H"
BCNG36	36L" x 2W" x 8H"
CUSTOM SIZES AVAILABLE	

