



HOSPITAL PRODUCTS



Imaging

Patient Room

Fall Prevention

Operating Room

Specialty Seating

Labor & Delivery

Emergency Room

Mattress & Bed Systems



Blue Chip

Medical Products Inc.

Tel: 845-369-7535 • Toll Free: 800-795-6115 • Fax: 845 369-7633

Website: www.bluechipmedical.com

"Setting the Standard"

OR & ER STRETCHER PAD & IMAGING TABLE PAD CONFIGURATIONS

CALL TO ORDER - 800 795-6115



MODEL 5003
Shown - One Section 3"
**See Options Below

MODEL 5023-GU
Shown - Two Section 3"
**See Options Below

MODEL 5033-GU
Shown - Three Section 3"
**See Options Below

GEL-FLEX™ OVERLAY

MODEL BC1000

CHOOSE A CONFIGURATION

Blue Chip manufactures OR & ER stretcher & imaging table pads for every configuration and manufacturer.

INCLUDING:

AMSCO® • MAQUET® • MIDMARK®
SHAMPAIGNE® • SIEMENS®
SKYTRON® • GE® • PHILLIPS®

THICKNESS OPTIONS



ONE SECTION

MODEL
5001 (1")
5002 (2")
5003 (3")
5004 (4")

TWO SECTION

MODEL	MODEL	MODEL	MODEL	MODEL
5021 (1")	5021-GU (1")	5021-GUF (1")	5021-GUBF (1")	
5022 (2")	5022-GU (2")	5022-GUF (2")	5022-GUBF (2")	
5023 (3")	5023-GU (3")	5023-GUF (3")	5023-GUBF (3")	
5024 (4")	5024-GU (4")	5024-GUF (4")	5024-GUBF (4")	

THREE SECTION

MODEL	MODEL	MODEL	MODEL
5031 (1")	5031-GU (1")	5031-GUF (1")	5031-GUBF (1")
5032 (2")	5032-GU (2")	5032-GUF (2")	5032-GUBF (2")
5033 (3")	5033-GU (3")	5033-GUF (3")	5033-GUBF (3")
5034 (4")	5034-GU (4")	5034-GUF (4")	5034-GUBF (4")


OR & ER STRETCHER & IMAGING TABLE CONFIGURATION CHART

SECTIONS:	1 • 2 • 3 • 4
THICKNESS:	1" • 2" • 3" • 4"
STOCK SIZES:	ALL MANUFACTURERS AND CONFIGURATIONS
CUSTOM SIZES:	CUSTOM SIZES AVAILABLE, PLEASE SEE ORDER FORM
CUT-OUTS:	PERINEAL CUT OUTS - CAN BE CUSTOMIZED GU (STANDARD CUT OUT BODY SECTION), GUF (CUT OUT FOOT SECTION), GUBF (CUT OUT BODY AND FOOT SECTION)
CORE OPTIONS:	FOAM (F) • VISCOTEC™ (VT) • SOFT-FLO™ GEL (SFG) • GEL-FLEX™ (GF) • GEL-FLEX™ CHANNEL (GFC) HONEYCOMB (SCHC) • SOFT-CUBE™ GEL (SCG) • GEL/FOAM - VISCO GEL/SWIRL (VGS) • VISCO GEL/BEADS (VGB) • AIR (GRZ)
COVER OPTIONS:	MEDICAL GRADE VINYL • CONDUCTIVE UL® • NON-CONDUCTIVE • VYVEX-I™ • VYVEX-II™ • VYVEX-III™ • NON-SKID
CORNER OPTIONS:	SQUARE • ROUND • CUT-OFF


CORE, COVER, CORNER & CLOSURE OPTIONS

CALL TO ORDER - 800 795-6115

CHOOSE A CORE




GEL-FLEX™
Flexible semi-solid Gel excellent for long surgical procedures. Provides superior pressure redistribution, available in .25", .50", .625", 1"



SOFT-FLO™ GEL
Aqueous Gel forms to encapsulate the bony prominences. Provides excellent pressure redistribution




GEL IMPREGNATED FOAM
Foam impregnated with gel provides excellent pressure redistribution. Encapsulates bony prominences



GEOMETRIC GEL
Hollow columns of solid gel in geometric shapes create pressure redistributing cells that deliver proper support. Special Soft-Cube™ heel gel for maximum protection against the incidence of heel ulcers


Blue Chip can customize OR, ER, Stretcher & Imaging Table pads for your application.



PPF™ - PRESSURIZED POLYMER FOAM
PPF™ is a polyurethane foam which greatly reduces the sag factor (deflection) enhancing the product's life & extending the uncompromising pressure redistributing qualities of our foam




VISCOTEC™ MEMORY FOAM
Heat sensitive visco elastic memory foam provides superior comfort and pressure redistribution qualities



GEL-FLEX™ CHANNEL
Combination of PPF™ foam and Gel-Flex™ provides excellent comfort and support for bariatric patients

GRZ™ NON-POWERED AIR

VYVEX-I COVER



GRZ™ NON-POWERED SELF ADJUSTING AIR

Nine self adjusting air bladders encapsulated in a high resiliency foam shell. 25" minimum width by any length, 4.95" height. Converts to powered system with optional pump

9 SELF ADJUSTING AIR BLADDERS

FOAM SHELL

OPTIONAL PUMP MODEL 4610



CHOOSE A COVER - CONDUCTIVE UL® OR NON-CONDUCTIVE

MEDICAL GRADE VINYL

Offered in a range of colors, Blue Chip's medical grade vinyl covers are fluid proof, anti-microbial and hypo-allergenic

NON-CONDUCTIVE VINYL

Non-conductive infection control covers are available in navy and black leatherette. Fluid-proof, anti-microbial and hypo-allergenic

CONDUCTIVE UL®

For surgical applications, covers are fluid-proof, anti-microbial & hypo-allergenic. Available in Vyvex-III™ stretch & vinyl

VYVEX I™

Nylon PVC-backed fabric, fluid proof, low shear, low friction, anti-microbial, stain resistant & washable

VYVEX II™


Nylon Butyl-backed fabric, fluid proof, low shear, low friction, anti-microbial, stain resistant and washable

VYVEX III™

4-way stretch urethane coated, fluid proof, low shear, low friction, anti-microbial, stain resistant and washable

NON-SKID BOTTOM

Durable non-skid pebbled bottom cover material can be added to help prevent accidental falls due to sliding or shifting during transfers



ALL COVERS ARE ANTI-MICROBIAL AND MEET CALIFORNIA TECHNICAL BULLETIN #117 FOR FIRE RETARDANCY

CHOOSE CORNER, CLOSURE & OPTIONS


CORNER STYLES



CLOSURE STYLES



OPTIONS



HEAD & NECK POSITIONERS

CALL TO ORDER - 800 795-6115



ADULT / PEDIATRIC / NEONATAL HORSESHOE

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck. Easy access for tubes

BC-1016	2" x 3" ID x 8" OD	Adult
BC-1016B	3" x 3" ID x 8" OD	Bariatric
BC-1019	1.25" x 2.25" ID x 5.50" OD	Ped
BC-1048	0.75" x 1.50" ID x 3.50" OD	Neo



ADULT / PEDIATRIC / NEONATAL DONUT

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck


BC-1005	2" x 3" ID x 8" OD	Adult
BC-1005B	3" x 3" ID x 8" OD	Bariatric
BC-1018	1.25" x 2.25" ID x 5.50" OD	Ped
BC-1047	0.75" x 1.50" ID x 3.50" OD	Neo



RICHARD TYPE HEAD REST

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck

BC-1036	3" x 9" x 8"
BC-1036-0	4.25" x 9" x 8" "T"



CONTOURED HEAD REST

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck


BC-1100	7" x 7" x 3.75"
---------	-----------------



PRONE SUPINE HEAD REST

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck


BC-1037	6" x 9" x 11" Cutout for Aen tubes
---------	------------------------------------



SUPINE & LATERAL HEAD PAD WITH DISH

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck

BC-1044	2" x 9" x 10" (Supine adult)
BC-1042	1" x 9" x 10" (Lateral adult)
BC-1055	2" x 6" x 7" (Supine pediatric)
BC-1056	1" x 6" x 7" (Lateral pediatric)
BC-1057	2" x 4" x 5" (Supine neo-natal)
BC-1058	1" x 4" x 5" (Lateral neo-natal)



OPHTHALMIC HEAD REST

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck for ophthalmic procedures


BC-1114	12" x 10" x 3.875"
---------	--------------------



ADULT HEAD CRADLE

DESC: PPF™ Blue Foam
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck. Designed for use in lateral supine or prone position


BC-11210	8" x 9" x 4"
----------	--------------



PEDIATRIC HEAD CRADLE

DESC: PPF™ Blue Foam
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck. Designed for use in lateral supine or prone position

BC-11213	6" x 7" x 3.5"
----------	----------------



SLOTTED HEAD CRADLE

DESC: PPF™ Blue Foam
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck. Designed to accommodate an endotracheal tube. Designed for use in lateral supine or prone position

BC-11212	8" x 8.5" x 4"
----------	----------------



ADULT & PEDIATRIC RING POSITIONERS

DESC: PPF™ Blue Foam
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck

BC-11202	1.5" x 7"
BC-11201	1.75" x 9"



ADULT & PEDIATRIC RING POSITIONERS

DESC: PPF™ Foam, Black conductive vinyl, flame retardant, water resistant seams, cleans easily
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck

BC-IR-10C-4ID	2" x 10" x 4" ID
BC-IR-9C-4ID	2" x 9" x 4" ID



ADULT & PEDIATRIC RING POSITIONERS

DESC: Charcoal foam, uncovered
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck

BC-IR-10-4ID-NC	2" x 10" x 4" ID
BC-IR-9-4ID-NC	2" x 9" x 4" ID



ADULT & PEDIATRIC RING POSITIONERS

DESC: PPF™ Foam, Black conductive vinyl, flame retardant, water resistant seams, cleans easily
 USE: Support, elevation, pressure redistribution and proper alignment of the head and neck

BC-IR-10C-6ID	2" x 10" x 6" ID
---------------	------------------


BODY POSITIONERS

CALL TO ORDER - 800 795-6115

FOAM BODY POSITIONERS

Cover: Grey Chemsafe™ Vinyl, Leatherette (See Color Chart)


Core: High Density Urethane foam



KNEE ELEVATOR

Shown with grey vinyl cover


KW-100	17L" x 10H" x 8W"
--------	-------------------



BODY ALIGNER

Shown with black vinyl cover

BA-100-CHAR	23W" x 8H" x 9.5D"
BA-100-V	23W" x 8H" x 9.5D"



POSITIONING BOLSTER ROLL

Shown with black vinyl cover

PBR-106V	10" x 6"
PBR-175V	17" x 5"
PBR-206V	20" x 6"
PBR-246V	24" x 6"

BED WEDGES


Shown with royal blue vinyl cover

Cover: Grey Chemsafe™ vinyl, Non-Conductive or Leatherette (See Color Chart), Conductive UL™ Material or Uncovered
Core: High Density PPF™ Polyurethane Foam

BW-475V	24"D x 26"W x 7.5"H
BW-410V	24"D x 26"W x 10"H
BW-412V	24"D x 26"W x 12"H

Custom sizes available

VINYL COLORS



Optional Non-skid Bottom available on all Body Positioners & Bed Wedges

GEL-FLEX™ BODY POSITIONERS

DESC: Gel-Flex™, durable, wipes clean, anti-bacterial non-allergenic

USE: Pressure redistribution, Positioning



PEDIATRIC CHEST ROLL

BC-1222	12" x 2" x 2"
BC-1233	12" x 3" x 3"
BC-0022	10" x 2" x 2"
BC-0033	10" x 3" x 3"
BC-0822	8" x 2" x 2"
BC-0833	8" x 3" x 3"
BC-0622	6" x 2" x 2"
BC-0633	6" x 3" x 3"
BC-0422	4" x 2" x 2"
BC-0433	4" x 3" x 3"

CHEST POSITIONER LEFT & RIGHT

BC-1045S	3" x 8" x 23/28"
----------	------------------

CHEST POSITIONER LEFT

BC-1045L	3" x 8" x 23/28"
----------	------------------



CHEST POSITIONER RIGHT

BC-1045R	3" x 8" x 23/28"
----------	------------------



ADULT CONTOURED CHEST ROLL

BC-1111	20" x 6.25" x 5"
BC-1112	16" x 6.25" x 5"
BC-1113	12" x 6.25" x 5"



ADULT ROUND CHEST ROLL

BC-1030	12" x 4" x 4"
BC-1031	12" x 6" x 6"
BC-1032	16" x 4" x 4"
BC-1033	16" x 6" x 6"
BC-1034	20" x 4" x 4"
BC-1035	20" x 6" x 6"



GRADUATED THORACIC POSITIONER

BC-1061	14" x 9" x 2.5"
---------	-----------------

ANKLE, KNEE, LEG, STIRRUP & SPECIALTY POSITIONERS


CALL TO ORDER - 800 795-6115

ABDUCTION PILLOWS

DESC: PPF™ Blue Foam, velcro straps

USE: Provides abduction of legs following hip surgery

BC-11323	12" x 18"	Small
BC-11322	17" x 22.5" H	Medium
BC-11321	18" x 25" H	Large



HEEL PROTECTORS

DESC: PPF™ Blue Foam, Velcro straps

USE: Heel protection

BC-11306	8" x 7"	Pair/Folded
----------	---------	-------------

HEEL CUP

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Maximum pressure redistribution by removing pressure from heels

BC-1046	1.75" x 3.25" x 6"
---------	--------------------




BC-1068	1.5" x 6" x 4.5"	Ped
BC-1069	1.5" x 8" x 4.5"	Ped
BC-1071	2.5" x 6" x 4.5"	Ped
BC-1075	1.5" x 12" x 4.5"	Adult
BC-1076	2" x 4.5" x 12"	Adult
BC-1039	1.5" x 12" x 16"	Adult
BC-1029	2.5" x 5" x 14"	Adult

SANDBAGS - ADULT AND PEDIATRIC

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Optimizes position and accessibility for surgical procedures




LATERAL POSITIONER

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Optimizes position and accessibility for surgical procedures

BC-1066	.50" x 20" x 28"
---------	------------------



WILSON FRAME PAD

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Prevents pressure and shearing during surgical procedures

BC-1043	.50" x 10" x 29"
BC-1043S	.50" x 10" x 29" Set



ANDREWS FRAME KNEE PAD

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Prevents pressure and shearing during surgical procedures

BC-1060	1" x 16" x 22"
---------	----------------

KNEE CRUTCH PAD

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Patients in lithotomy position, urology or gynecology procedures

BC-1025	.375" x 9" x 11"	Small
BC-1009	.375" x 11" x 17"	Medium
BC-1008	.50" x 20" x 20"	Large
BC-1038	.50" x 11" x 27"	Extra large

BARIATRIC KNEE CRUTCH

BC-1052	1" x 9" x 11"	Small
BC-1053	1" x 11" x 17"	Medium
BC-1054	1" x 20" x 20"	Large



BARIATRIC LAP BELT PAD

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Lateral and seating during surgical/medical procedures

BC-1050	.25" x 6" x 36"	Medium
BC-1051	.25" x 7" x 48"	Large



ANKLE STRAP STIRRUP PAD

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic. Comes in pairs.

USE: Patients in lithotomy position, urology or gynecology procedures

BC-1082	.375" x 1.875" x 12.5"	Small
BC-1083	.375" x 1.875" x 17.375"	Large



CANDY CANE STIRRUP PAD

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Patients in lithotomy position, urology or gynecology procedures

BC-1027	2" I.D. X 25" Tube
---------	--------------------

ARM & ELBOW POSITIONERS

CALL TO ORDER - 800 795-6115

ARMBOARD

DESC: PPF™ Foam, black non-conductive vinyl cover

USE: Supportive padding during procedures

OPT: Thickness: 2", 3", 4"

Core: PPF™ Foam

Core: PPF™ Foam, Viscotec™ Memory Foam & Gel-Flex™

Core: PPF™ Foam, Soft-Flo™ & Viscotec™ Memory Foam


Core: PPF™ Foam & Viscotec™ Memory Foam

Cover: Black conductive vinyl, Vyvex-I™, Vyvex-II™, Vyvex-III™


Please specify cover, thickness or core options when ordering

BCAB-624	6" x 24" (2", 3", 4")
BCAB-626	6" x 26" (2", 3", 4")



CUTAWAY
VISCOTEC™
SOFT-FLO™
PPF™



CUTAWAY
GEL-FLEX™
VISCOTEC™
PPF™



OPTIONAL CORE THICKNESS



ARM CRADLE

DESC: PPF™ Blue Foam

USE: Supportive padding during surgery
Contoured design keeps the arm in place while the notches allow the device to bend with the patient's arm

BC-1331	5" x 23" x 3"
---------	---------------



CONTOURED ARMBOARD

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Positioning & pressure relief during surgical procedures

BC-1084	2.25" x 5.75" x 20" Medium
BC-1088	2.25" x 5.75" x 16" Small




CONVOLUTED ARM PAD

DESC: Convoluted, PPF™ Blue Foam

USE: Supports and provides protection

BCPF-220	18" x 9" x 2"
BCPF-221	20" x 6" x 2"
BCPF-900	24" x 10" x 2"



ARMBOARDS

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Positioning & pressure relief during surgical procedures

BC-1017	.25" x 4" x 15" Small
BC-1118	.25" x 5" x 20" Medium
BC-1006	.25" x 5" x 25" Large

BARIATRIC ARMBOARDS

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Positioning & pressure redistribution for surgical procedures

BC-1007	.375" x 6" x 26"- XLarge
BC-1049	.375" x 6" x 15"- XLarge Short



ELBOW PAD

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Positioning & pressure distribution for surgical procedures

BC-7710	3.25" OD x 1.50" x .75"
---------	-------------------------



SQUARE ULNAR PAD

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Provides protection of the ulnar nerve & elbow

BC-1020	3.75" x 11" x 17" Square
---------	--------------------------



OVAL ULNAR PAD

DESC: Gel-Flex™ Durable, wipes clean, anti-bacterial, non-allergenic

USE: Provides protection of the ulnar nerve & elbow

BC-1021	3.75" x 6" x 18" Oval
---------	-----------------------



ELBOW PAD PROTECTOR

DESC: PPF™ Blue Foam, velcro straps, non-sterile, radiolucent

USE: Provides protection of the ulnar nerve & elbow

BC-1301	6" x 15.5" x .5" Pair
---------	-----------------------

HOSPITAL & PATIENT ROOM SUPPORT SURFACES

CALL TO ORDER - 800 795-6115

LABOR & DELIVERY, PEDIATRIC

L & D MATTRESS



DESC: PPF™ Polyurethane high density foam with choice of anti-microbial cover
 USE: Birthing Bed

BC-BIR	SPECIFY BRAND NAME & MODEL #	CUSTOM SIZES AVAILABLE
--------	------------------------------	------------------------

HOSPITAL, INSTITUTION, PRISON & DRUG REHAB



INSTITUTIONAL MATTRESS

DESC: PPF™ Polyurethane high density foam, Medical Grade Vinyl, Sealed seams, Inverted lock-stitch closure. Standard cover grey Chemsafe™ vinyl.
 USE: For Long Term and Acute Care Facilities, Psychiatric, Institutional & Prison

BC-PSY-3680	36" x 80" x 6"	CUSTOM SIZES AVAILABLE
-------------	----------------	------------------------

SPECIALTY SEATING, DIALYSIS, CHEMO, DENTAL



BASSINET MATTRESS

DESC: Fire-retardant PPF™ high-compression polyurethane foam, covered in Grey Chemsafe™ Vinyl, Sealed seams
 OPT: Thickness - 1", 2", Cover - Vinyl, Vyvex-I™, Nylon

BC-BASS	SPECIFY SIZE	CUSTOM SIZES AVAILABLE
---------	--------------	------------------------



WARMER MATTRESS

DESC: PPF™ Polyurethane foam with choice of Grey Chemsafe™ Vinyl or Vyvex-I Nylon cover, Sealed seams
 USE: Neo-natal warming mattress

BC-WARM	SPECIFY SIZE	CUSTOM SIZES AVAILABLE
---------	--------------	------------------------



PEDIATRIC CRIB MATTRESS

DESC: Fire-retardant high-compression PPF™ polyurethane foam, covered in Grey Chemsafe™ Vinyl, Sealed seams
 USE: Hospital and in Home
 OPT: Thickness - 1", 2", 3", Cover - Vinyl, Vyvex-I™ Nylon

BC-PED	26" x 50" x 4"	CUSTOM SIZES AVAILABLE
--------	----------------	------------------------



GERI GEL RECLINER OVERLAY

DESC: PPF™ Polyurethane foam, Gel, Vyvex-III™ cover
 Non-skid bottom, Supports 250 lbs
 USE: Dialysis & Chemo recliner overlay

6200NS	18" x 70" x 1.75"
--------	-------------------



GERI GEL BARIATRIC RECLINER OVERLAY

DESC: Gel-Flex™ Channel, Vyvex-III™ Cover / Non-skid bottom
 Supports 800 lbs
 USE: Bariatric recliner overlay

BCGG-6000-BAR	29" x 70" x 4"	CUSTOM SIZES
---------------	----------------	--------------



COMFORT SUPREME DENTAL CHAIR OVERLAY

DESC: MEDFLEX VPF™, ViscoTec™ Memory foam, Vyvex-III™ Cover / Non-skid bottom,
 USE: Comfort & Pressure redistribution for dental procedures

6600	Navy	64" x 22" - 14" x 3"
6600B	Black	64" x 22" - 14" x 3"
6600R	Rose	64" x 22" - 14" x 3"
6600T	Teal	64" x 22" - 14" x 3"



GERI CHAIR ALTERNATING PAD & AC/DC PUMP

DESC: Alternating Air Bladders, Vyvex-III™ cover
 USE: Geriatric Recliner comfort & pressure redistribution
 *Chair not included

9700GR	21" x 70" x 3.5"
--------	------------------

PRESSURE REDISTRIBUTION POWERED & NON-POWERED MATTRESS & BED SYSTEMS

CALL TO ORDER - 800 795-6115



AIRUS™
Air Fluidized Therapy
Bed System

Model A210

Model A210	32" x 84"	PDAC E0194
------------	-----------	------------

- First Group 3 Portable Air Fluidized Therapy, AFT Bed
- Full Length Air Fluidized Therapy
- Sterilized Support Beads provides clean, dry therapeutic surface
- Hand Held Control Unit
- Medical Grade Breathable Nylon top cover
- Cover is anti-microbial and fire retardant
- Foam Headrest
- Whisper Quiet Pump
- IV Pole & Triangle Pull Bar
- Heavy Duty Trolley for ease of portability



POWER-PRO™
Alternating Pressure
with Low Air Loss Mattress System

Model 9500

Model 9500	36" x 80" x 9"	PDAC E0277
------------	----------------	------------

- For the prevention & treatment of stages I - IV pressure ulcers
- Supports 650 lbs.
- Bariatric Models to support up to 1,000 lbs.
- 18 - 9" deep ventilated bladders
- Available widths 36", 42" & 48"
- Vyvex III™ Urethane coated multi-stretch Low shear, fluid proof top cover
- Cover is anti-microbial and fire retardant




POWER-TURN™
Lateral Rotation Mattress System

Model 9800

Model 9800	36" x 80" x 9"	PDAC E0277
------------	----------------	------------

- For the prevention & treatment of stages I - IV pressure ulcers
- Full Turn, Left Turn, Right Turn, 25° Small Turn & 40° Big Turn
- Supports 650 lbs.
- Bariatric Models to support up to 1,000 lbs.
- 18 - 9" deep ventilated bladders
- Available widths 36", 42" & 48"
- Vyvex III™ Urethane coated multi-stretch Low shear, fluid proof top cover
- Cover is anti-microbial and fire retardant



GRZ™
Self Adjusting Air Mattress

Model GRZ 3676




Optional Pump
Model 4610

Dynamic Air Bladders

Model GRZ 3676	36" x 76" x 7"	PDAC E0277
----------------	----------------	------------

- For the prevention & treatment of stages I - IV pressure ulcers
- Supports 550 lbs.
- Bariatric Models to support up to 1,000 lbs.
- 9 Individual Dynamic Self Adjusting Air bladders
- High Resilience Foam Shell
- Available widths 36", 42", 48", 54", 60"
- Available lengths 76", 80", 84"
- Vyvex I™ Low shear nylon, fluid proof top cover
- Cover is anti-microbial and fire retardant
- Optional Alternating Pressure Pump



THE PRE-VENT II SERIES®
Therapeutic pressure redistributing foam mattresses

Shown without Cover


Model BCEL-

BCEL	Specify Size
------	--------------

- For the prevention & treatment of stages I - II pressure ulcers
- Standard 6" supports 350 lbs.
- Bariatric Models to support up to 1,000 lbs.
- Tri-Laminate Design
- Available widths 36", 39", 42", 48" & 54"
- Available lengths 76", 80", 84"
- Heel Wrap™ PPF™ foam for the prevention of heel ulcers
- Vyvex II™ Low shear nylon, fluid proof top cover
- Foams & cover materials are anti-microbial and fire retardant

FALL PREVENTION & PATIENT SAFETY PRODUCTS

CALL TO ORDER - 800 795-6115



SIDE RAIL BUMPER PAD ENVELOPE STYLE

DESC: Fire-retardant high-compression PPF™ polyurethane foam, covered in anti microbial vinyl. Protects top rail

USE: Protection from injury due to banging or slipping through a bed rail

BCM-211-ENV	11" x 70" x 1"
BCM-211-2-ENV	11" x 70" x 2"
BCM-215-ENV	15" x 60" x 1"
BCM-215-2-ENV	15" x 60" x 2"
BCM-216-ENV	11" x 30" x 2"



SIDE RAIL BUMPER PAD WITH SEE THROUGH WINDOW

DESC: Fire-retardant high-compression PPF™ polyurethane foam, covered in anti microbial vinyl

USE: Protection from injury due to banging or slipping through a bed rail

BCM-211-W	11" x 70" x 1"
BCM-211-2-W	11" x 70" x 2"
BCM-215-W	15" x 60" x 1"
BCM-215-2-W	15" x 60" x 2"
BCM-216-W	11" x 30" x 2"



DEFENDER EDGE FALL MAT

DESC: Fire-retardant closed cell foam, Grey Chemsafe™ vinyl cover / non-skid bottom.

USE: Fall Protection, edge allow wheelchair access

BCEDG-2566	25" x 66" x 1.5"
BCEDG-3072	30" x 72" x 1.5"




FALL MAT

DESC: Fire-retardant high-compression PPF™ polyurethane foam, medical grade vinyl cover

Available in 2", 3" & 4" in a bi-fold, tri-fold or non-fold design. Non skid bottom

USE: Fall Protection

BCFM-2460	24" x 60" x 2"
BCFM-2472	24" x 72" x 2"
BCFM-3072	30" x 72" x 2"
BCFM-3672	36" x 72" x 2"



SIDE RAIL BUMPER PAD WITH VELCRO STRAPS

DESC: Fire-retardant high-compression PPF™ polyurethane foam, covered in anti microbial vinyl

USE: Protection from injury due to banging or slipping through a bed rail

BCM-211	11" x 70" x 1"
BCM-211-2	11" x 70" x 2"
BCM-215	15" x 60" x 1"
BCM-215-2	15" x 60" x 2"
BCM-216	11" x 30" x 2"




NO-GAP SIDE GUARD

DESC: Fire-retardant high-compression PPF™ polyurethane foam covered in anti microbial vinyl

USE: Protect patients from accidental slippage through gaps between the mattress and side rail

BCNG-72	72" x 2" x 8"
BCNG-36	36" x 2" x 8"



MATTRESS EXTENDER

DESC: Fire-retardant high-compression PPF™ polyurethane foam covered in anti microbial vinyl

USE: Used to extend mattress for the needs of a taller patient. Can be customized to fit any bed

BCME-126	6" x 36" x 6"
----------	---------------



WEDGE GUARD

DESC: Fire-retardant high-compression PPF™ polyurethane foam, covered in anti microbial vinyl

USE: Side rail gap protector & patient positioner
Patients are safe from falls and accidental slippage through gaps

BWG-36	HALF RAIL 36"
BWG-70	FULL RAIL 70"

CHOOSE COLORS



Grey
(Standard)




Navy
(Optional)



Green
(Optional)



Blue
(Optional)



Tan
(Optional)

CUSTOM MEDICAL CART, EQUIPMENT COVERS & PILLOWS

CALL TO ORDER - 800 795-6115



Shown with optional routing pocket



Shown with optional clear front flaps

Blue Chip manufactures customized cart covers for all purposes & equipment manufacturers

Utility Carts
Crash Carts
Linen Carts
Medical Carts
Case Carts
Isolation Carts

Guard Rail Carts
Multi-Shelf Carts
Emergency Carts
Anesthesia Carts
Narcotics Cabinets
Specialty Carts

Medical Surgical Carts
Wire Shelf Utility Carts
Multi-Purpose Carts
Medical Laptop Carts
Healthcare Carts
and More...

Optional: Blue Chip will custom imprint your cart cover upon request. Custom imprint can be silk screened or embroidered.


HOSPITAL PATIENT PILLOWS

DESC: Hypo-allergenic, high lofted, fire retardant fiber fill. Available covers: vinyl, breathable vinyl, navy nylon and poly-cotton. The nylon and vinyl covers are anti-microbial, and meet California Technical Bulletin #117 for fire-retardancy.



2026V	20" x 24"	Fire Retardant with Anti-Microbial Vinyl Cover
2026B	20" x 24"	Fire Retardant with Breathable Anti-Microbial Vinyl Cover
2026N	20" x 24"	Fire Retardant with Anti-Microbial Nylon Cover

CHOOSE CLOSURE



VELCRO - STANDARD




ZIPPER - OPTIONAL

CHOOSE OPTIONS

ROUTING POCKET

Specify Left Side, Right Side or Front Center




CHOOSE 400 DENIER NYLON COLOR



CHOOSE CART COVER FABRIC


BALLISTIC NYLON

1000 Denier - Heavy Duty. Available in Royal



VINYL

See Vinyl Color Chart below




400 DENIER NYLON

Available in a large variety of colors. See chart



CHOOSE VINYL COLORS



MANUFACTURING • WAREHOUSE • SALES • SERVICE

Blue Chip has a complete manufacturing facility. Our products are manufactured to medical & military standards and can be customized for your specific application.



INVENTORY

A huge inventory of raw materials and stocked manufactured products means the products you need are in-stock and ready for quick shipment from our New York factory



EXTENSIVE
MATERIAL
CHOICE &
CUTTING
CAPABILITY



SALES, SERVICE & SUPPORT

Blue Chip provides the industries finest sales, service & technical support. All products are quality tested to ensure effectiveness and reliability before they ship to you