

Air-Pro® Series

Alternating Pressure Pad & Pump Overlay Mattress Systems

AIR-PRO® SERIES - OVERLAY MATTRESS SYSTEM

The Air-Pro® Plus & Elite Alternating Pump and Overlay Mattress Systems include a foldaway bed hanger bracket. This feature keeps the pump off the floor and reduces any accidental damage. Our easy access filter located on the bottom of the pump reduces any chance of dust particles from entering the pump compressor motor. We include a heavy duty, deluxe bubble style pad with head and foot tuck under flaps. The adjustability of the Elite allows this system to be used for heavier patients. Air-Pro® meets California Technical Bulletin #117 for fire retardancy.



Model 4400 system shown

AIR-PRO® ALTERNATING PRESSURE PAD AND PUMP OVERLAY MATTRESS SYSTEMS

- Air-Pro® Plus **Model 4200**
- Air-Pro® Elite **Model 4400**

 **Blue Chip**
Medical Products, Inc.
"Setting the Standard"

AIR-PRO®

SPECIFICATIONS

AIR-PRO® SERIES - OVERLAY MATTRESS PAD

Material: Medical grade PVC with Anti-microbial, fire-retardant treatments.

Design: Bubble Style with more touch points for effective pressure redistribution.

Model 4810 shown



	Air-Pro® Pad Solo	Air-Pro® Pad Deluxe	Air-Pro® Pad Laser	Air-Pro® Pad Multi
Model:	4800	4810	4820	4830
Size:	35" x 79" x 2.5"	35" x 79" x 2.5"	35" x 79" x 2.5"	35" x 79" x 3.75"
End Flaps:	No	Yes	Yes	Yes
Ventilated:	No	No	36 Lasered Holes	36 Lasered Holes
Pump:	4201/4401	4201/4401	4401S	4401S
Heavy Duty:	No	Yes	Yes	Yes

MODEL 4200 MODEL 4201

System - 4201 Pump w/ 4810 Overlay Pad Pump Only

Dimensions: 9L" x 4W" x 3.5H"
Weight: 4 lbs.
Electrical: 110V 60 hz
Listings:
Cycle Time: 10 Minutes
Air Output: 5.6 Liters
Power Cord: 12 ft.

*All models available in 220-240V; 50-60Hz



MODEL 4400 MODEL 4401 MODEL 4401S

System - 4401 Pump w/ 4810 Overlay Pad Pump Only Pump Only (Models 4820 & 4830)

Dimensions: 10L" x 5W" x 4H"
Weight: 5 lbs.
Electrical: 110V 60 hz
Listings:
Cycle Time: 10 Minutes
Air Output: 6.5 Liters
Air Output (4401S): 8 Liters
Power Cord: 12 ft.

*All models available in 220-240V; 50-60Hz



WARRANTY

Pump: 1 year limited
Mattress: 6 months limited

PDAC Approved
Medicare E0180



Blue Chip
Medical Products, Inc.

7-11 Suffern Place, Suffern NY 10901

Tel: (845) 369-7535 • Fax: (845) 369-7633 • Toll Free: (800) 795-6115
www.bluechipmedical.com

"Setting the Standard"

©Blue Chip Medical Products, Inc. 2009

Alternating Pressure Mattresses
 Alternating Pressure Pads & Pumps
 Low Air-Loss Mattress Systems
 Lateral Rotation Mattress Systems
 Gel Mattresses
 Gel Mattress Overlays
 Static-Air Mattress Overlays
 Water Mattress Overlays
 Wheelchair Cushions
 Patient Lift Slings
 Side Rail Bumpers
 Cart Covers
 OR Replacement Table Pads
 Replacement Stretcher Pads
 Cubicle Curtains
 Hospital Pillows