



Therapeutic Wheelchair Back System

Patients who spend long periods of time in their wheelchairs need more back support and comfort. The Amara® Therapeutic Wheelchair Back System is the solution. Blue Chip has designed the perfect complement to the Amara® component based therapeutic wheelchair seating system. This creates a therapeutic seating environment customized to the individual patient.



Amara® Therapeutic Wheelchair Back System improves the comfort and support of your existing wheelchair back by eliminating the hammock effect. Amara provides proper patient positioning to maintain spinal and shoulder alignment and superior pressure redistribution to prevent skin breakdown.

Amara® Therapeutic Wheelchair Back System comes with either a fixed or adjustable lumbar support and is easy to install. In addition, Amara® is available in a Dual-Zone design which combines the unmatched support of Medflex® foam with the therapeutic pressure redistribution qualities of Viscotec® foam.

The Amara® line of wheelchair back cushions (E2611-E2612) is designed for a patient who has a wheelchair and meets the Medicare coverage criteria. The use of these products provide better sitting tolerance and comfort for patients requiring a wheelchair.



"Setting the Standard"

WHEELCHAIR BACKS

AMARA® - FIXED BACK

- High density Medflex® foam base construction
- Supports the back to provide a comfortable upright sitting position
- Rigid base with a fixed lumbar support creates an ultra-lightweight back with a low profile design
- Vyvex® III non-shear incontinent top cover
- Non-Skid bottom cover
- Cover is washable, anti-microbial and fluid resistant
- Adjustable safety straps properly attach to chair
- Back is fire retardant and meets Cal. Tech #117

MODEL	SIZE	COVER
AMR-BK-2-1618	16" x 18"	Vyvex® III Top / Non-Skid Bottom
AMR-BK-2-1818	18" x 18"	Vyvex® III Top / Non-Skid Bottom
AMR-BK-2-2020	20" x 20"	Vyvex® III Top / Non-Skid Bottom
AMR-BK-2-2220	22" x 20"	Vyvex® III Top / Non-Skid Bottom
AMR-BK-2-2420	24" x 20"	Vyvex® III Top / Non-Skid Bottom

Supports 250 lbs.

CUSTOM AND BARIATRIC SIZES AVAILABLE, BASED ON PATIENT'S WEIGHT

Limited 18 Month WARRANTY



AMARA® - DUAL-ZONE VISCOTEC® BACK

- High density Medflex® foam base construction with therapeutic top layer of Viscotec® for added comfort and pressure redistribution
- Supports the back to provide a comfortable upright sitting position
- Vyvex® III non-shear incontinent top cover
- Non-Skid bottom cover
- Cover is washable, anti-microbial and fluid resistant
- Adjustable safety straps properly attach to chair
- Back is fire retardant and meets Cal. Tech #117



MODEL	SIZE	COVER
AMR-BK-1-1618	16" x 18"	Vyvex® III Top / Non-Skid Bottom
AMR-BK-1-1818	18" x 18"	Vyvex® III Top / Non-Skid Bottom
AMR-BK-1-2020	20" x 20"	Vyvex® III Top / Non-Skid Bottom
AMR-BK-1-2220	22" x 20"	Vyvex® III Top / Non-Skid Bottom
AMR-BK-1-2420	24" x 20"	Vyvex® III Top / Non-Skid Bottom

Supports 250 lbs.

CUSTOM AND BARIATRIC SIZES AVAILABLE, BASED ON PATIENT'S WEIGHT

Limited 18 Month WARRANTY

AMARA® - ADJUSTABLE BACK

- High density Medflex® foam base construction
- Supports the back to provide a comfortable upright sitting position
- Rigid base with an adjustable lumbar support creates an ultra-lightweight back with a low profile design.
- Vyvex® III non-shear incontinent top cover
- Non-Skid bottom cover
- Cover is washable, anti-microbial and fluid resistant
- Adjustable safety straps properly attach to chair
- Back is fire retardant and meets Cal. Tech #117



MODEL	SIZE	COVER
AMR-BK-3-1618	16" x 18"	Vyvex® III Top / Non-Skid Bottom
AMR-BK-3-1818	18" x 18"	Vyvex® III Top / Non-Skid Bottom
AMR-BK-3-2020	20" x 20"	Vyvex® III Top / Non-Skid Bottom
AMR-BK-3-2220	22" x 20"	Vyvex® III Top / Non-Skid Bottom
AMR-BK-3-2420	24" x 20"	Vyvex® III Top / Non-Skid Bottom

Supports 250 lbs.

CUSTOM AND BARIATRIC SIZES AVAILABLE, BASED ON PATIENT'S WEIGHT

Limited 18 Month WARRANTY



7-11 Suffern Place, Suffern NY 10901

Tel: (845) 369-7535 • Fax: (845) 369-7633 • Toll Free: (800) 795-6115

www.bluechipmedical.com

'Setting the Standard'