


# Air-Pro™ Series

*Alternating Pressure Pad & Pump Overlay Mattress Systems*

## AIR-PRO™ SERIES - OVERLAY MATTRESS SYSTEM

The Air-Pro™ Plus & Elite Alternating Pump and Overlay Mattress Systems include a foldaway bed hanger bracket. This feature keeps the pump off the floor and reduces any accidental damage. Our easy access filter located on the bottom of the pump reduces any chance of dust particles from entering the pump compressor motor. We include a heavy duty, deluxe bubble style pad with head and foot tuck under flaps. The adjustability of the Elite allows this system to be used for heavier patients. Air-Pro™ meets California Technical Bulletin #117 for fire retardancy.

Shown with  
4810 pad



Model 4400 system shown

- |                             |                           |
|-----------------------------|---------------------------|
| • Air-Pro™ <sup>Plus</sup>  | Model 4200 with flaps     |
| • Air-Pro™ <sup>Plus</sup>  | Model 4200S without flaps |
| • Air-Pro™ <sup>Elite</sup> | Model 4400 with flaps     |
| • Air-Pro™ <sup>Elite</sup> | Model 4400S without flaps |



# AIR-PRO™

## SPECIFICATIONS

### AIR-PRO™ SERIES - OVERLAY MATTRESS PAD

**Material:** Medical grade PVC with Anti-microbial, fire-retardant treatments.

**Design:** Bubble Style with more touch points for effective pressure redistribution.



Model 4800 Pad  
without flaps



Model 4810 Pad  
with flaps

	Air-Pro™ Pad Solo	Air-Pro™ Pad Deluxe	Air-Pro™ Pad Laser	Air-Pro™ Pad Multi
<b>Model:</b>	4800	4810	4820	4830
<b>Size:</b>	35" x 79" x 2.5"	35" x 79" x 2.5"	35" x 79" x 2.5"	35" x 79" x 3.75"
<b>End Flaps:</b>	No	Yes	Yes	Yes
<b>Ventilated:</b>	No	No	36 Lasered Holes	36 Lasered Holes
<b>Pump:</b>	4201/4401	4201/4401	4401SS	4401SS
<b>Heavy Duty:</b>	Yes	Yes	Yes	Yes
<b>Patient Weight:</b>	250	250	250	350

#### MODEL 4200

#### MODEL 4200S

#### MODEL 4201

**System - 4201 Pump w/ 4810 Overlay Pad**


**System - 4201 Pump w/ 4800 Overlay Pad**

**Pump Only**

**Dimensions:** 9L" x 4W" x 3.5H"

**Weight:** 4 lbs.

**Electrical:** 110V 60 hz

**Listings:** 

**Cycle Time:** 10 Minutes

**Air Output:** 5.6 Liters

\*All models available in 220-240V; 50-60Hz



Model 4201

#### MODEL 4400

#### MODEL 4400S

#### MODEL 4401

#### MODEL 4401SS

**System - 4401 Pump w/ 4810 Overlay Pad**

**System - 4401 Pump w/ 4800 Overlay Pad**

**Pump Only**

**Pump Only (to be used with Laser & Multi pad only)**

**Dimensions:** 10L" x 5W" x 4H"

**Weight:** 5 lbs.

**Electrical:** 110V 60 hz

**Listings:** 

**Cycle Time:** 10 Minutes

**Air Output:** 6.5 Liters

**Air Output (4401SS):** 8 Liters

\*All models available in 220-240V; 50-60Hz



Model 4401



Model 4401SS  
To be used with  
Laser & Multi pad only

PDAC Approved  
Medicare E0181

#### WARRANTY

**Pump:** 1 year limited  
**Mattress:** 6 months limited



Blue Chip

Medical Products, Inc.

7-11 Suffern Place, Suffern NY 10901

Tel: (845) 369-7535 • Fax: (845) 369-7633 • Toll Free: (800) 795-6115

www.bluechipmedical.com

"Setting the Standard"