


Support Surfaces

Pressure Redistribution Mattresses



Alternating Pressure Systems
with Low Air Loss



Lateral Rotation Systems
with Low Air Loss



Convertible Mattress Systems



Homecare & Institutional Mattresses



Alternating Pressure Pads & Pumps




 **Blue Chip**
Medical Products, Inc.

7-11 Suffern Place, Suffern NY 10901
Toll Free: (800) 795-6115 • Fax: (845) 369-7633
www.bluechipmedical.com

"Setting the Standard"

ALTERNATING PRESSURE WITH LOW AIR LOSS



Model 4600

APOLLO 3-PORT™

The Apollo 3-Port™ Alternating Pressure Low Air Loss Mattress System combines the most advanced technology and clinically effective features in a portable system. Provides pressure redistribution for the prevention and treatment of Stages I-IV pressure ulcers. Model 4600SRS with Safety Side Rail System.

Model:	4600
Supports:	350 lbs
Size:	35"W x 80"L x 8"H
Bladders:	18 - 8" air cells, ventilated, low air loss
Fire Code:	Cal. Tech Bulletin #117
Cover:	***Vyvex III™
Warranty:	1 Year Limited
PDAC:	E0277



Model 9600



Ventilated
Foot Section
9600 & 4600

SUPREME AIR™

The Supreme Air™ Alternating Pressure Low Air Loss Mattress System combines the most advanced technology with clinically effective static and alternating features in a portable system. Provides pressure redistribution for the prevention and treatment of Stages I-IV pressure ulcers. Includes low pressure alert and audible power outage alarm. Model 9600SRS with Safety Side Rail System.

Model:	9600 (35"), 9642 (42")
Supports:	350 lbs
Size:	35"W x 80"L x 8"H, 42"W x 80"L x 8"H
Bladders:	18 - 8" air cells, ventilated, low air loss
Fire Code:	Cal. Tech Bulletin #117
Cover:	***Vyvex III™
Warranty:	1 Year Limited
PDAC:	E0277



Optional Digital Pump

POWER PRO™ ELITE

Combines two clinically effective therapies: Alternating Pressure with Low Air Loss in a portable blower based system. Provides proper pressure redistribution for the prevention and treatment of Stages I-IV pressure ulcers. The new Power-Pro VW Elite bariatric model fits 35", 39", 42" and 48" bed frames. An optional pulsating low air loss therapy feature is available. Model 9500SRS with Safety Side Rail System.



Model 9500

Model:	9500
Supports:	650 lbs. - 1,000 lbs.
Size:	36"W x 80"L x 9"H
	Available in 39"W, 42"W, 48"W
Bladders:	18 - 9" air cells, ventilated, low air loss
Fire Code:	Cal. Tech Bulletin #117
Cover:	***Vyvex III™
Warranty:	1 Year Limited
PDAC:	E0277

TRADEWIND™ SERIES

The Tradewind™ series offers therapists a choice of four dual mode pumps. All Tradewind™ pumps feature Alternating Pressure and Static modes as standard. The Alternating Pressure Mode provides the therapeutic benefit of reactive hyperemia and redistributes pressure well below 32mmHg, while the Static Mode provides gentle low air loss with constant pressure redistribution and support. This allows nursing to perform patient care.



TRW-1000
Mechanical pump

TRW-3000
Digital Pump

TRW-ATS
Digital Pump

TRW-ATS-RC
Digital Pump
Remote Control




Model:	TRW-1000, TRW-3000, TRW-ATS, TRW-ATS-RC
Supports:	350 lbs. - 850 lbs.
Size:	35"W x 80"L x 8"H avail. in 42"W, 48"W
Bladders:	17 - 8" air cells, ventilated, low air loss
Fire Code:	Cal. Tech Bulletin #117
Cover:	***Vyvex III™
Warranty:	1 Year Limited
PDAC:	E0277

***Vyvex III™ is multi-stretch knit urethane, low shear, fluid and stain proof, anti-microbial and fire retardant.

AIR SYSTEMS

GRZ™ - NON-POWERED CONVERTIBLE SYSTEM

The GRZ™ is a non-powered, pressure redistribution, self-adjusting, convertible mattress. GRZ™ helps to provide proper pressure redistribution to reduce interface pressures by dynamically adjusting to body weight and movement. The GRZ™ air cells are encased in a Medflex™ pressure redistribution PPF™ foam shell to redistribute pressure where it is needed most. Optional pump converts GRZ™ to alternating air therapy. Maintenance Free.



Model:	GRZ™
Supports:	up to 1,000 lbs based on size
Size:	36"W x 80"L x 7"H
	36"W x 76"L x 7"H also available in
	39"W, 42"W, 48"W, 54"W 60"W
Bladders:	9 self adjusting air cells
Fire Code:	16 CFR 1633
Cover:	Vyvox II™ Top, Medical grade vinyl bottom
Warranty:	5 years Limited
PDAC:	E0277 (requires alternating pump)

LATERAL ROTATION WITH LOW AIR LOSS

POWER TURN™ ELITE

The Power-Turn™ Elite Mattress System combines two clinically effective therapies: Lateral Rotation with Low Air Loss in a portable blower based system provides proper pressure redistribution. Lateral Rotation Therapy reduces the risk of pulmonary complications for patients who are immobile. A Full, Left or Right Turn mode is offered with either a 25° Small Turn or a 40° Big Turn.



Model:	9800
Supports:	650 lbs. - 1,000 lbs.
Size:	36"W x 80"L x 9"H avail. in 42"W, 48"W
Bladders:	18 - 9" air cells, ventilated, low air loss
Fire Code:	Cal. Tech Bulletin #117
Cover:	***Vyvox III™
Warranty:	1 Year Limited
PDAC:	E0277


MATTRESS OPTIONS

SAFETY SIDE RAIL SYSTEM - SRS

SRS Safety Side Rail System can be added to all Alternating Pressure Mattresses for additional patient safety. To reduce accidental falls.

OPTIONAL INFLATOR

Blower pump can be used through the CPR Valve to inflate the mattress in under four minutes.



RAPID-AIR™ and RAPID-AIR PLUS™


Alternating Air Floatation Mattress Systems with Gentle Low Air Loss provides pressure redistribution for the prevention and treatment of Stages I-IV pressure sore ulcers. The model 4300P Rapid-Air Plus™ features a PPF™ 2.5" pressure redistribution foam base that protects the patient from bottoming out during power outage. Pediatric model 4300PM. Model 4300SRS with Safety Side Rail System.




Model:	4300 / 4300P / 4300PM
Supports:	275 lbs.
Size:	35"W x 80"L x 5.5"H (standard)
	25"W x 45"L x 3.5"H (pediatric)
Bladders:	18, 5.5" air cells
Fire Code:	Cal. Tech Bulletin #117
Cover:	***Vyvox III™
Warranty:	1 Year Limited
PDAC:	E0277

***Vyvox III™ is multi-stretch knit urethane, low shear, fluid and stain proof, anti-microbial and fire retardant.


THERAPEUTIC FOAM SLEEP SYSTEMS




SUSPENSION
Model BCSS II



SUPREME
Model BCSP II



ELITE
Model BCEL II




PRO
Model BCPR II


PRE-VENT II SERIES®

The Pre-Vent II Series® features five pressure redistribution foam mattresses. The Pre-Vent II Series® incorporates Blue Chip Medical's exclusive Sphere™, Reflect™ SCT™ and Multi-Flex™ foam contour designs. These react to body weight displacement across the entire surface to prevent and treat pressure sores. Blue Chip's unique construction together with our Med-Flex™ PPF™ foams, react instinctively to a compression set, allowing for a fully automatic, totally zoned pressure-redistribution surface. All models feature the exclusive Heel Wrap™ PPF™ foam heel section. Custom sizes and bariatric mattresses to support patients up to 1000 lbs available.

Model:	Suspension, Supreme, Elite, Pro, RC
Supports:	350 lbs, Bariatric 650 lbs - 1000 lbs. Elite only
Size:	36"W x 80"L x 6"H, (76"L & 84"L)
Addl. Widths	39"W, 42"W, 48"W, 54" & 60"
Core:	Med-Flex™ PPF™ - Reflect™ SCT™
Cover:	BCSS - Vyvex III™, BCSP - Vyvex II™
	BCEL - Vyvex II™, BCPR - Vyvex I™, BCRC - Vyvex I™
Warranty:	BCSS 8 Yr., BCSP 7 Yr., non pro-rated
	BCEL 6 Yr., BCPR 6 Yr., BCRC 5 Yr., non pro-rated



RC II
Model BCRC II



Perimeter Safety Border
(PB) Standard



Shown with Optional
Safety Side Rail System (SRS)
(must be added at time of ordering)



PRE-VENT II
Shown In Cover

HEEL WRAP
Sphere



Reflect




PRE-VENT II ELITE™ BARIATRIC EXPANDABLE

Developed especially for the rental market, the Pre-Vent II Elite™ Expandable bolsters allow for use with 39", 42", 48" and 54" width bedframes and 84" & 88" lengths. The mattress is standard at 36" with 6" side bolster to expand to 42" & 48" and 4" head & foot bolsters to extend the length of the mattress to 84" & 88" Mattress is also available in 42" expandable to 54". The Vyvex-II® nylon butyl-back top cover incorporates sleeves that allow for easy access, insertion or removal of the bolsters. Gel-Flex™ can be added as an option.

Model:	BCEL-EXP
Supports:	Bariatric 650 lbs - 1000 lbs.
Size:	36"W x 80"L x 7"H, (84", 88")
Addl. Widths	39"W x 42"W x 48"W, 54"W & Custom Sizes
Core:	Med-Flex™ PPF™ - Reflect™ SCT™
Warranty:	6 Yr. non-pro-rated, Gel-Flex 7 yr non-pro-rated

42" expandable to 54". The Vyvex-II® nylon butyl-back top cover incorporates sleeves that allow for easy access, insertion or removal of the bolsters. Gel-Flex™ can be added as an option.



Model BCEL-EXP-GLFLX
Shown with Gel-Flex

*FLAMMABILITY

All Mattresses Meet Federal Flammability Standard 16 CFR 1633, California Technical Bulletin #117 & #129, 603, Life Safety Code 101 and City of Boston BF 1X-11.

THERAPEUTIC FOAM SLEEP SYSTEMS

RELIEF-CARE™ PRO

The Mattress System uses advanced Surface Cut Technology (SCT™) to create a clinically effective support surface. SCT™ enables Blue Chip to produce a therapeutic dual-zone PPF™ Pressurized Polymer Foam support surface with varying features and properties. Blue Chip's design absorbs and redistributes patient's weight to provide greater comfort, maximum support and proper pressure redistribution around bony prominences. Custom sizes available.



Model BCHTM

Model:	BCHTM
Supports:	250 lbs.
Size:	36"W x 80"L x 5.5"H (76"L, 84"L)
Core:	PPF™ Pressurized Polymer Foam
Cover:	Vyvox I™ Top, Medical grade vinyl bottom
Fire Code:	16 CFR 1633
Warranty:	1 Year Limited
PDAC:	E0184

HOME-CARE ELITE™

Designed for use in patient homes, hospitals, and institutions. The inner core is a high density polyurethane foam that is covered with an anti-microbial vinyl material and closed with a self-ventilating zipper. *Mattress meets Federal Flammability Standards. Custom sizes available. 1 Year Limited Warranty.



Model BCHM


Model:	BCHM
Supports:	250 lbs.
Size:	36"W x 80"L x 5.5"H
Core:	High Density Polyurethane Foam
Cover:	Chemsafe™ Vinyl

INSTITUTION AND BEHAVIORAL MATTRESS

The inner core is a high density polyurethane foam that is covered with a heavy duty anti-microbial vinyl cover. The cover is sealed with an inverted lock-stitch to protect against infiltration using a super-strong nylon thread. *Mattress meets Federal Flammability Standards. Available in grey, navy or Clear Safe® vinyl. We also offer an optional mesh vent and/or built-in pillow. Custom sizes available. 1 Year Limited Warranty.



Model BCPR II-3680-BM



Model BCPR II-3680-BMBIP
Shown with built in pillow

Model:	BCPR II-3680-BM
Supports:	250 lbs.
Size:	36"W x 80"L x 6"H
Core:	High Density Polyurethane Foam
Cover:	Heavy Duty Chemsafe™ Vinyl

CRIB & BASSINET MATTRESS

Designed for use in hospitals, patient homes and institutions. The inner core is a high density polyurethane foam that is covered with Blue Chip's Chemsafe™ anti-microbial vinyl cover and closed with a self ventilating zipper or can be ordered with sealed seams. *Mattress meets Federal Flammability Standards. Custom sizes available. 1 Year Limited Warranty.



Model BC-PED



Model BC-BASS

Model:	BC-PED
Thickness:	4"
Size:	Custom
Core:	High Density Polyurethane Foam
Cover:	Chemsafe™ Vinyl

Model:	BC-BASS
Thickness:	1', 1.5", 2"
Size:	Custom
Core:	High Density Polyurethane Foam
Cover:	Chemsafe™ Vinyl

OEM REPLACEMENT MATTRESS COVERS & REPLACEMENT COVERS WITH SAFETY SIDE RAIL SYSTEM

Ideal for refurbishing a worn mattress without incurring the expense of replacement. The OEM replacement cover is made from Chemsafe™ medical grade vinyl or optional Vyvox-I™, Vyvox-II™ or Vyvox-III™ top covers. They are fluid proof and offer a low-shear, low-friction surface and effectively protect against a wide spectrum of micro-organisms, including MRSA, mold, mildew, staining and odors. These covers feature Blue Chip's exclusive Barrier-Lok™ infection control system. Custom sizes are available.



Model:	BCOEM976
Size:	36"W x 76"L x 6"H
Model:	BCOEM980
Size:	36"W x 80"L x 6"H

Model:	BC-CVR-SRS-EZ-3680
Size:	36"W x 80"L x 6"H

*Vyvox I™ is a soft nylon with PVC back, **Vyvox II™ is 70 denier nylon, with butyl back, ***Vyvox III™ is a multi-stretch knit urethane. All materials are low-shear fluid and stain resistant, anti-microbial and fire retardant.

The SRS-EZ replacement mattress cover with raised foam safety side rail improves patient safety of any foam mattress. Top cover is made with Vyvox-I™ nylon and a medical grade vinyl bottom. Complies with FDA Patient Entrapment Guidelines. Optional in Vyvox-II™ or Vyvox-III™ multi-stretch.



SRS-EZ

*FLAMMABILITY

All Mattresses Meet Federal Flammability Standard 16 CFR 1633, California Technical Bulletin #117 & #129, 603, Life Safety Code 101 and City of Boston BF 1X-11.

MATTRESS OVERLAYS

GEL-PRO™ MATTRESS OVERLAY & GEL-PRO™ MATTRESS

The Gel-Pro™ mattress overlay Model 6103 offers maximum protection to guard against the incidence of pressure sores. Our viscous gel, migrates through sectioned bladders to provide proper immersion and envelopment of the bony prominences with over 80% surface coverage. A high density foam base supports the entire mattress overlay for maximum weight distribution and durability. Available in bariatric, full and queen sizes.

The Gel-Pro line of mattresses Model 6305, is ideal for the Nursing Home and Homecare Markets. The 3 section, 5.5 inch mattress utilizes a high density, high compression foam base designed to extend the life of the mattress and maximize its clinical benefit. The mattress is constructed using the same materials as the Gel Overlay but with heavy duty use in mind. The zippered cover allows for sections to be removed for easy portability and storage. Materials are fluid and stain resistant, anti-microbial, and fire-retardant. Also available in bariatric, full and queen sizes.



Model 6103

Model:	6103
Supports:	250 lbs.
Size:	35"W x 76"L x 3"H
Bariatric:	42", 48", 54" & 60" x 80" X 4" (weight based)
Cover:	Vyvox I™ Top, Medical Grade vinyl bottom
Core:	Viscous Gel / Foam
Fire Code:	Cal. Tech Bulletin #117
Warranty:	1 Year Limited
PDAC:	E0185



Model 6305

Model:	6305
Supports:	300 lbs
Size:	36"W x 80"L x 5.5"H
	42"W x 80"L x 6"H, 48"W x 80"L x 6"H
Core:	Gel / High Density Foam
Cover:	Vyvox I™ Top, Medical grade vinyl bottom
Warranty:	1 Year Limited



STAT-H2O™, STAT-GEL™

The Stat-H2O™ water mattress helps to redistribute pressure while providing a cooler more comfortable sleeping surface. The internal baffle design balances the patient in the center of the mattress while eliminating the traditional water mattress "rocking" effect. Stat-H2O™ easily converts to STAT-GEL™, a full gel mattress overlay with the addition of our special polymer gel powder.



Model 5500

Model:	STAT-H2O™ 5500, STAT-GEL™ 5700
Supports:	250 lbs.
Size:	32"W x 72"L x 3.5"H
Bladders:	Medical Grade Vinyl with Baffles
Warranty:	6 Months
PDAC:	5500 E0198, 5700 E0185

STAT-AIR™

The mattress is constructed using a pliable medical grade vinyl that contours to the patient's body. Heat build up and moisture are greatly reduced, due to the air holes throughout the mattress. For patient safety adjustable anchor straps secure Stat-Air™ to an existing mattress. Quick release valve for rapid deflation during medical emergencies. The mattress is anti-microbial and fire retardant.



Model 5300

Model:	5300
Supports:	250 lbs.
Size:	35"W x 72"L x 3.5"H
Mattress:	Medical Grade Vinyl
Warranty:	6 Months
PDAC:	E0197



Electric Inflator
Model 9623



Manual Inflator
Model 5301

AIR-PRO™ PLUS, AIR-PRO™ ELITE

The Air-Pro™ Plus Model 4200 & Elite Model 4400 Alternating Pump and Overlay Mattress Systems feature a heavy duty deluxe bubble style pad with head and foot tuck under flaps. Pads are available with or without end flaps and lasered air holes for ventilation. The adjustability of the Elite allows this system to be used for heavier patients. Multipad with high output 4400S pump supports patient weights to 350 lbs.



Model 4200
Model 4400

Model:	Plus 4200, Elite 4400
Pads:	4800, 4810, 4820, 4830
Supports:	250 lbs., 4830 supports 350 lbs.
Size:	35"W x 79"L x 2.5"H
Mattress:	Bubble style with flaps
Fire Code:	Cal. Tech Bulletin #117
Warranty:	Pump 1 Year, Mattress 6 months
PDAC:	E0181