

Gel Pro® Wheelchair Cushions

GEL-PRO® ELITE - SKIN PROTECTION GEL WHEELCHAIR CUSHIONS

Our gel bladder contains an aqueous, highly viscous gel which provides for pressure redistribution, proper immersion and envelopment around the bony prominences. The bladder is surrounded by a high density foam shell that is covered using Vyvex-I®, a fluid-proof, low shear nylon top with a medical grade vinyl bottom. All cover materials are anti-microbial and meet California Technical Bulletin #117 for fire retardancy. An optional solid seat insert is available.



PDAC Approved
 Medicare E2603 - less than 22"W
 Medicare E2604 - 22"W or Greater

Cover:

Top: Vyvex I® - Low shear, fluid proof and anti-microbial nylon cover
Bottom: Medical grade vinyl

Core:

Top: Med-Flex® VPF® Foam
Center: Ultra soft viscous gel bladder
Base: Med-Flex® VPF® high density foam

MODEL	SIZE	COVER	THICKNESS
6204NA	16" x 16"	Vyvex-I® Top / Vinyl Bottom	3"
6204	17" x 17"	Vyvex-I® Top / Vinyl Bottom	3"
6204A	18" x 16"	Vyvex-I® Top / Vinyl Bottom	3"
6204FA	20" x 16"	Vyvex-I® Top / Vinyl Bottom	3"
6204VFA	22" x 16"	Vyvex-I® Top / Vinyl Bottom	3"
6204VVFA	24" x 16"	Vyvex-I® Top / Vinyl Bottom	3"
6204C	18" x 18"	Vyvex-I® Top / Vinyl Bottom	3"
6204CC	20" x 18"	Vyvex-I® Top / Vinyl Bottom	3"
6204CCFA	22" x 18"	Vyvex-I® Top / Vinyl Bottom	3"
6204CCVFA	24" x 18"	Vyvex-I® Top / Vinyl Bottom	3"
6204D	20" x 20"	Vyvex-I® Top / Vinyl Bottom	3"

Supports 250 lbs.

Custom and Bariatric sizes available, based on patient's weight

Limited 18 Month Warranty

WHEELCHAIR CUSHIONS

GEL-PRO® LOW PROFILE - LOW PROFILE GEL WHEELCHAIR CUSHION

An excellent choice where a low profile cushion is necessary. Our segment gel bladder contains an aqueous, highly viscous gel which provides for pressure redistribution, proper immersion and envelopment around the bony prominences. The Gel Pro Low Profile is covered using Vyvex-1®, a fluid-proof, low shear nylon top with a medical grade vinyl bottom. All cover materials are anti-microbial and meet California Technical Bulletin #117 for fire retardancy. An optional solid seat insert is available.

MODEL	SIZE	COVER	THICKNESS
6201NA	16" x 16"	Vyvex-1® Top / Vinyl Bottom	1.75"
6201	17" x 17"	Vyvex-1® Top / Vinyl Bottom	1.75"
6201A	18" x 16"	Vyvex-1® Top / Vinyl Bottom	1.75"
6201FA	20" x 16"	Vyvex-1® Top / Vinyl Bottom	1.75"
6201VFA	22" x 16"	Vyvex-1® Top / Vinyl Bottom	1.75"
6201VFVA	24" x 16"	Vyvex-1® Top / Vinyl Bottom	1.75"
6201C	18" x 18"	Vyvex-1® Top / Vinyl Bottom	1.75"
6201CC	20" x 18"	Vyvex-1® Top / Vinyl Bottom	1.75"
6201CCFA	22" x 18"	Vyvex-1® Top / Vinyl Bottom	1.75"
6201CCVFA	24" x 18"	Vyvex-1® Top / Vinyl Bottom	1.75"
6201D	20" x 20"	Vyvex-1® Top / Vinyl Bottom	1.75"
Supports 250 lbs.			

Custom and Bariatric sizes available, based on patients weight
 Limited 12 Month Warranty



Cover:

Top: Vyvex 1® - Low shear, fluid proof and anti-microbial nylon cover

Core:

Top: Med-Flex® VPF® High Density Foam
Center: Ultra soft viscous gel bladder
Base: Medical Grade Vinyl

GEL PRO® DLX - CONVOLUTED GEL WHEELCHAIR CUSHION

This versatile reversible gel cushion utilizes the same gel formulation and bladder design as used in the Gel-Pro Elite cushion. The low shear nylon cover is anti-microbial and fluid-proof to protect the inner core. The Gel Pro DLX Profile is covered using Vyvex-1®, a fluid-proof, low shear nylon top with a medical grade vinyl bottom. All cover materials are anti-microbial and meet California Technical Bulletin #117 for fire retardancy.

Zippered Cover:

Cover: Vyvex 1® - Low shear, fluid proof and anti-microbial nylon cover

Core:

Top: Convoluted Med-Flex® VPF® Foam
Center: Ultra soft viscous gel bladder
Base: Convoluted Med-Flex® VPF® Foam



MODEL	SIZE	COVER	THICKNESS
606NA	16" x 16"	Vyvex-1®	3"
606	17" x 17"	Vyvex-1®	3"
606A	18" x 16"	Vyvex-1®	3"
606FA	20" x 16"	Vyvex-1®	3"
606VFA	22" x 16"	Vyvex-1®	3"
606VFVA	24" x 16"	Vyvex-1®	3"
606C	18" x 18"	Vyvex-1®	3"
606CC	20" x 18"	Vyvex-1®	3"
606CCFA	22" x 18"	Vyvex-1®	3"
606CCVFA	24" x 18"	Vyvex-1®	3"
606D	20" x 20"	Vyvex-1®	3"
Supports 250 lbs.			

Custom and Bariatric sizes available, based on patients weight
 Limited 12 Month Warranty



7-11 Suffern Place, Suffern NY 10901
 Tel: (845) 369-7535 • Fax: (845) 369-7633 • Toll Free: (800) 795-6115
www.bluechipmedical.com

"Setting the Standard"