

Home-Care^{Elite}®

Homecare & Institutional High-Density Foam Mattress

Designed for use in Nursing Homes, Hospitals, and Institutions. The inner core is a high density, fire-retardant polyurethane foam that is covered with an anti-microbial, fire-retardant material meeting Federal Flammability Standard 16 CFR 1633. The mattress is closed with a self-ventilating zipper, and can be special ordered for use in prisons and psychiatric areas utilizing an inverted lock-stitch closure.

Single layer of high-density foam encased in a zippered medical grade vinyl cover providing superior patient comfort.



The Home-Care Elite® is covered with a zippered medical grade vinyl material. The cover is fluid proof, anti-microbial, fire-retardant and stain resistant. The cover material and foam core meet Federal Flammability Standard 16 CFR 1633 for fire retardancy. Mattress can be manufactured to meet California Technical Bulletin #129 and Boston Fire Code 1X-11 when specified.



Durability, resiliency and patient comfort make the Home-Care Elite® the perfect mattress for homecare and institutional applications.

MODEL	SIZE	THICKNESS
BCHM 760	36" x 76"	5.5"
BCHM 780	36" x 80"	5.5"
BCHM 784	36" x 84"	5.5"
BCHM 742	42" x 80"	7"

CUSTOM AND BARIATRIC SIZES AVAILABLE

WARRANTY: 12 MONTHS



Blue Chip

Medical Products, Inc.

7-11 Suffern Place, Suffern NY 10901

Tel: (845) 369-7535 • Fax: (845) 369-7633 • Toll Free: (800) 795-6115

www.bluechipmedical.com

"Setting the Standard"